

Metzeler Mattress

Instructions for Use



With the purchase of this Metzeler Mattress, you have decided on a brand product from Metzeler Schaum GmbH “made in Germany”.

For over 60 years, Metzeler have been producing high-quality mattresses, and places great emphasis on supporting your individual requirements for a restful night's sleep.

In order to enjoy your product for as long as possible, please observe the following points:

In new mattresses, what is known as “new product odour” can occasionally occur, but generally dissipates within just a few days. However, in some cases it can take up to eight weeks for it to disappear completely. All foam types used have been tested to the Okotex 100 Product Group 1 standard.

As mattresses are made of soft material, it is possible that slight alterations of length and width may occur because of transport and storage. A deviation of up to 20mm in length and width is possible. Should greater deviations be found, please contact your specialist dealer before using the mattress.

Metzeler mattresses adapt ideally to your body. For this reason, it is possible that, during the initial phase of use, a slight hollow may form, particularly in the hip area, caused by your body weight, in particular, this depression can be evident after sleep. The performance characteristics are not affected as a result. Mattresses with an integrated seating edge are firmer in the outer zones than in the middle.

Mattress Cleaning

We recommend airing the mattress regularly so that moisture absorbed during the night can be released into the room. Pay attention to ensuring the bedroom is well ventilated with an ideal humidity level (40–60%).

Your Metzeler mattress has been tested to withstand standard cleaning products in line with RKI, VAH, and IHO lists. Metzeler recommends following cleaning instructions as directed by cleaning product manufacturers procedures.

The cover is machine washable. Remove sticky label from washable covers before laundering. Wash covers separately. Use mild detergent and the gentle cycle on your machine. Please note care instructions in the care label. Not suitable for tumble dryer.

Mattress Maintenance

Double-sided mattresses should ideally be turned every time bed linen is changed, by both flipping the mattress over and swapping the head and foot ends around to ensure the mattress is used as evenly as possible. The mattress has been designed for use in a horizontal position. Permanent use with extreme bed articulation can damage the mattress.

We recommend the mattress is inspected for any damage or significant wear during every linen change and replace the mattress or cover if deemed necessary.